



ExoBand

MOVEO S.R.L.

USE AND MAINTENANCE MANUAL

Hip orthosis

ExoBand

Version: 2.2 - Date of latest version of the manual May 2022

Table of Contents

DECLARATION OF CONFORMITY	2
MANUFACTURER'S INFORMATION.....	2
AUTHORIZED TECHNICAL SERVICE	2
THE MANUAL.....	2
Conventions.....	3
SYMBOLS.....	4
WARRANTY	4
GENERAL SAFETY RULES	5
Warning.....	5
DESCRIPTION OF THE MEDICAL DEVICE	5
TECHNICAL DATA.....	6
INTENDED USE AND UNINTENDED USE OF THE DEVICE	7
Sphere of application and intended use.....	7
Unintended use.....	7
Limitations of the medical device	7
RESIDUAL RISKS.....	7
Adverse effects	8
Contraindications.....	8
HANDLING	8
USE.....	9
Storage	9
Application.....	9
Precautions.....	11
MAINTENANCE	14
Safety	14
Periodic maintenance	14
Washing.....	16
Special maintenance	16
Disposal	16

DECLARATION OF CONFORMITY

The medical device described in this documentation is sold with the declaration of conformity in compliance with the regulations in force in Europe.



NOTE

BEFORE USING THE DEVICE, MAKE SURE THAT THE DECLARATION OF CONFORMITY HAS BEEN DELIVERED TOGETHER WITH IT.

MANUFACTURER'S INFORMATION

Company name	Moveo S.r.l.
Registered office:	Via Monsignor Fortin 37/38, 35128 Padua - Italy
VAT no.	05236760285
Tel. Office	+39 049 261 44 27
Mobile	+39 3914590627
E-mail	info@moveowalks.com
Website	www.moveowalks.com

AUTHORISED TECHNICAL SERVICE

Technical operations on the hip orthosis may be carried out exclusively by authorized and qualified personnel of Moveo S.r.l. or its representative.

THE MANUAL



NOTE

READ THIS MANUAL CAREFULLY BEFORE CARRYING OUT ANY OPERATION WITH THE DEVICE.

This manual contains instructions for the use and maintenance of the ExoBand Hip Orthosis. The manual is composed of sections, each of which addresses different topics, subdivided into chapters and paragraphs.

The table of contents lists all the topics covered in this manual. Each page is numbered progressively. This manual contains all the instruction for the proper use, maintenance, and storage of the hip orthosis during its lifespan.

This manual contains confidential information and must not be provided to third parties, even partially, for any use and in any form, without the prior written consent of the manufacturer.

Moveo S.r.l. declares that the information contained in this manual is consistent with the technical and safety specifications of the medical device to which the manual refers. A certified copy of this manual is contained in the device's technical file kept by Moveo S.r.l.

Moveo S.r.l. does not recognize any documentation that has not been produced, released, or distributed by Moveo S.r.l. itself or by one of its authorized representatives.

This manual, like the entire technical file, will be kept by the manufacturer for the period provided for by law.

During this period, customers may request a copy of the documentation provided with the product at the time of purchase.

The entire technical file will be available to supervisory authorities, which may request a copy thereof, for the period provided for by law.

After this period, whoever uses and/or handles the product must make sure that both the product and the documentation comply with the laws in force.

Conventions

In order to facilitate comprehension of the topics covered, the manual utilizes graphic and typographic symbols and conventions as described below.

Graphic warning symbols



NOTE

NOTES CONTAIN CRITICAL INFORMATION HIGHLIGHTED APART FROM THE TEXT TO WHICH THEY REFER



WARNING







WARNINGS INDICATE PROCEDURES WHICH MUST BE STRICTLY FOLLOWED, AS TOTAL OR PARTIAL NON-COMPLIANCE MAY CAUSE DAMAGE TO THE MEDICAL DEVICE AND MAY EXPOSE THE PATIENT TO DANGER.



DANGER

INDICATIONS OF DANGER REGARD PROCEDURES WHICH MUST BE FOLLOWED WITHOUT EXCEPTIONS, AS TOTAL OR PARTIAL NON-COMPLIANCE MAY CAUSE DAMAGE OR INJURY TO THE PATIENT OR OTHERS NEARBY.

SYMBOLS

	Lot number
	Date of manufacture
	Manufacturer
	Agent
	Read the manual before use
	CE mark

WARRANTY

The terms of the warranty, listed in full in the purchase contract, are only valid if the medical device is used as intended and exclusively for the purposes for which it was designed.

Except for the routine and special maintenance operations described in the **MAINTENANCE** section and carried out in accordance with the indicated procedures, any repair or modification made to the device by the user or by unauthorized companies will void the warranty.

The warranty does not extend to damage caused by incompetence or negligence in the use of the device, or as a result of poor or lacking maintenance.

The products we sell are covered by a warranty under the following conditions:

1	<i>The warranty is valid for a period of twelve (12) months for juridical persons and twenty-four (24) months for natural persons, starting from the date indicated on the invoice. The warranty on parts subject to wear (mechanical components) is valid for a period of six (6) months from the date indicated on the invoice.</i>
2	<i>The manufacturer undertakes to replace all defective parts at its own discretion, only after inspecting the device/the parts in question.</i>
3	<i>Shipping costs are always borne by the purchaser if the terms of the warranty are not met.</i>
4	<i>During the warranty period, the replaced components become the property of the manufacturer.</i>
5	<i>Only the original purchaser who has complied with the maintenance instructions contained in this manual can benefit from this warranty. Our warranty liability expires when the original purchaser relinquishes ownership of the device, or if the device is modified.</i>
6	<i>The warranty does not cover damage resulting from excessive stress (use of the device after a fault or defect has been found, use of unsuitable operating methods, as well as failure to observe the use and maintenance instructions).</i>
7	<i>The manufacturer is not liable for difficulties arising in the resale or use of the device abroad as a result of the provisions of the laws in force in the country where the device was sold.</i>

Note: If it is deemed necessary to use the warranty, please indicate the following information:

1	Type
2	Date of purchase (submit the purchase receipt)
3	Detailed description of the problem

GENERAL SAFETY RULES

Warning

This section describes the general safety rules to be observed during any operation performed with the medical device. The operating procedures described in the following sections must be carried out in accordance with the operating methods indicated and the general safety provisions provided in this section.



NOTE

THE MANUFACTURER SHALL NOT BE LIABLE FOR ACCIDENTS OR DAMAGE CAUSED BY INAPPROPRIATE USE OF THE PRODUCT, OR BY EVEN PARTIAL FAILURE TO COMPLY WITH THE SAFETY RULES AND OPERATING PROCEDURES DESCRIBED IN THIS MANUAL.

Failure to comply with the instructions for the use and maintenance of the medical device contained in this manual also voids the warranty.

DESCRIPTION OF THE MEDICAL DEVICE

WARNINGS

The medical device must be prescribed and used under medical supervision. The orthopedic specialist is the expert of reference both for the application and for information relating to safe use. To guarantee the efficacy, tolerability and proper operation of the product, it must be applied and adjusted with the utmost accuracy. Any modification of the structure of the device or adjustment of the latter must be prescribed by a doctor and implemented by an orthopedic specialist. In hypersensitive individuals, direct contact with the skin may cause redness or irritation. If you experience pain, swelling, or any other abnormal reaction, contact your doctor immediately.



The medical device is designed exclusively to support human walking. Any unintended use is not recommended (such as using it when travelling in a car, riding a motorbike or bicycle, or even in the water, etc.). Any damage resulting from improper use is the responsibility of the customer.

The device consists of a belt, two leg loops and two tensioners which together with a toothed guide join all these elements. The device is available in 5 sizes (XS, S, M, L, and XL).

Class I Medical device.

TECHNICAL DATA

Belt and mechanism:

- Airshell 100% Polyester
- AL
- FE 10B21
- Cordura® 40% Cotton Lycra and 60% Polyamide
- Nylon 618 fiber
- Polyvinyl Chloride (PVC) 82% - Polyethersulfone (PES) 18%
- Spandex 10% - Polyester 90%
- Brass with nickel barrel plating finish
- Polyester 70% - Spandex 30%
- Polyoxymethylene (POM)
- 100% polypropylene
- Polyamide/polyurethane resin
- 87% Polyamide - 13% acrylic
- Nylon 37% - Polyester 33% - Rubber and Silicone 30%
- NY6 Buckles (Polimid B AV Natural HF)
- Adhesive interlining 100% Polyamide
- Tough Resin FLTOTL05 / PA12
- Tough RPU 70

Leg loops:

- Airshell 100% Polyester
- Cordura® 40% Cotton Lycra and 60% Polyamide
- Acrylic (LOCTITE 243)
- Galvanized steel
- Aluminum alloy sleeves
- Polyamide/polyurethane resin
- 87% Polyamide - 13% acrylic
- Nylon 37% - Polyester 33% - Rubber and Silicone 30%
- Aluminum (Nickel) and plastic
- Polyester 75% - Rubber 25%
- Polyvinyl Chloride (PVC) 82% - Polyethersulfone (PES) 18%
- Stainless steel
- Adhesive interlining 100% Polyamide
- Polyester 70% - Spandex 30%
- Tough Resin FLTOTL05 / PA12
- Tough RPU 70
- Polytetrafluoroethylene (PTFE)
- Dyneema® SK78 - HT Polyester

Features: Belts designed to adapt the product to different human morphologies.

INTENDED USE AND UNINTENDED USE OF THE DEVICE

Sphere of application and intended use

The medical device must be used exclusively as:

- a walking aid for the elderly and people with motor impairments caused by disabling diseases and/or health conditions
- an aid to improve improper or pathological posture.



NOTE

THE MANUFACTURER IS NOT LIABLE FOR ACCIDENTS OR DAMAGE CONSEQUENTIAL TO THE UNINTENDED USE OF THE PRODUCT. ANY UNINTENDED USE OF THE DEVICE WILL ALSO VOID THE WARRANTY.

Unintended use

The medical device is not designed for any use other than those described in the paragraph "**SPHERE OF APPLICATION AND INTENDED USE**". Furthermore, the following uses are absolutely forbidden:

- use of the medical device by individuals who have not read the manual
- use by children
- use by pregnant women
- use in direct contact with the skin.

A medical prescription is recommended for use by individuals suffering from serious diseases such as musculoskeletal system tumors, neuromuscular pathologies, spinal pathologies, hernias, etc.

Limitations of the medical device

- The medical device is not able to cure serious spinal deformations (dysmorphisms) or other dysmorphisms (it can only reduce the outcomes).
- The medical device can be used for long periods of time, that is, until it shows signs of wear or damage that are such as to affect its structure, safety, and operation. In this case and when in doubt, consult the seller immediately.
- The medical device can reach its maximum efficiency after a gradual period of use according to the psychophysical conditions of the user.

RESIDUAL RISKS

During the design phase Moveo S.r.l. conducted an in-depth risk analysis of the device. This analysis demonstrated that some risks cannot be eliminated due to their nature. These residual risks are addressed in this manual and indications regarding how to prevent them are provided throughout the same. It is therefore important that anyone who uses or handles the device has read and understood the manual beforehand.



NOTE

THE MANUFACTURER IS NOT LIABLE FOR ACCIDENTS OR DAMAGE CONSEQUENTIAL TO ANY UNINTENDED OR NEGLIGENT USE OF THE PRODUCT.

In particular:

- Do not make any modification to the medical device. Any damage resulting from the use of the device if this has been improperly modified by an unauthorized person releases the manufacturer from any liability.
- Keep this manual in good condition and easily available as this is necessary for the proper and safe use of the hip orthosis.
- If the structure of the device has sharp edges following an accidental collision or fall or a part worn enough to make it hazardous, contact an authorized assistance service and follow the instructions provided. +++

Adverse effects

In the analysis and experimentation phase, no adverse effects associated with the use of the medical device were detected. Sores, redness, numbness or tingling can occur in case of improper use of the medical device.

Contraindications

- The length of time the medical device is used must be proportional to the psychophysical conditions of the user.
- The medical device may cause muscle pain and/or fatigue when used without observing the indications contained in "Limitations of the medical device."
- Elastic bands, as such, if stretched and then released can cause pain or injury. They must therefore be carefully applied and handled to prevent their uncontrolled release (as with any elastic material).
- Keep to and ensure compliance with the instructions provided above in section "Unintended use of the medical device".

HANDLING

Carefully check the integrity of the device and its parts upon receipt. If there are any signs of damage, deformation, impact, or missing parts, inform the manufacturer before proceeding with any subsequent operation. The medical device must always be handled with care to avoid damaging it and making it unusable and/or dangerous. It can be handled manually without any problem.

Procedure for reporting damage or defects

If any damage or defects are found, stop the assembly process and report the nature of the damage or defect found to Moveo S.r.l.'s customer support team.


USE

Storage

- Store the hip orthosis indoors, away from the elements and humidity.
- Store the hip orthosis away from heat sources, open flames, and direct sunlight.
- Store the leg loop tensioners properly.
- Store the product in a box and/or drawer.

Use

- 1 Position the belt at the height of the anterior iliac crests (ASIS)
- 2 Carefully tighten the Velcro strap of the belt keeping the sliding buckles (A) at the height of the ASIS
- 3 Close the side release buckles (B) and tighten the fabric straps (C), tightening the belt and anchoring it to the pelvis without creating excessive compression
- 4 Once the belt is tightened, check that the sliding buckles (A) are aligned over the ASIS
- 5 Align the rear end of the belt with the spinal column
- 6-7 Once the sliding buckles (A) are in place, fasten the strip with the pressure buckle (D) located on the side of the belt
- 8 Hook the tensioner (E) of the leg loops to the toothed guide (F) of the belt. Place the leg loops close to the kneecap.
- 9 Fasten the leg loops starting from the two inner fabric wings (G) and tighten them carefully. Fasten the leg loops on the soft part of the Velcro strap (H).
- 10 Pull and attach the elastic covering strip (K) to the soft part of the Velcro strap (H).
- 11-12 Adjust the tension of the mechanism (L) by lifting the metal tab of the tensioner (E) by 90 degrees until you hear it "click."





Warning:


- The medical device must be worn over clothing.
- Trousers with raised rivets and/or buttons can damage the belt.
- Do not overtighten the mechanism (L) as this may make the device uncomfortable.
- To loosen or release the mechanism (L), press the button in the center of the tensioner (E) and push it downwards as shown in the figures below.

**Precautions**

Improper use may damage the device.

Front and rear belt positioning

The belt must be positioned at the correct height as illustrated in section "USE" of this manual. The proper and improper positioning of the belt is shown below.



Excessive tension of the mechanism


Excessive tension of the mechanism results in the irregular alignment of the belt as shown in the figures below. If the device is in this configuration, decrease the tension of the mechanism.



Inserting the tensioner into the toothed guide


Insert the tensioner into the toothed guide by proceeding as shown in the figures, holding it on its sides.

Do not squeeze the spring tensioner as, once inserted, it could prevent the tensioner from sliding and damage the sliding mechanism.



Aligning the leg loops

Always make sure that the leg loops are properly aligned with the belt as shown in the figure. Misalignment may cause the cable to break and consequently make the device unusable.



MAINTENANCE




Safety

It is absolutely necessary to read the manual before carrying out any maintenance operation.


Periodic maintenance

The hip orthosis can be periodically cleaned to eliminate dirt and bad smells. To wash it, remove the mechanism (L) following the procedures below:

- 1- Remove the screw cap (M) located on the upper part of the leg loop (Figure 1 - 1A).
- 2- Gently turn the screw cap (M) counterclockwise as shown in Figure 2 until reaching position 2A.
- 3- Remove the elastic element and store it in a safe place, taking care not to damage the cable (Figure 3).
- 4- Situation after removing the elastic element.



- 5- After washing the device, screw the screw cap (M) back into its seat, making sure to insert it properly. Insert the screw cap (M) in the direction shown in Figure 5A.
- 6- Rotate it clockwise until reaching the position shown in Figure 6A.
- 7- Final result.



Washing

Machine wash in cold water on a delicate cycle after removing the elastic components from the leg loops as described above. Leave to dry on a drying rack, away from heat sources. Do not use dryer. Do not use fabric softener. Do not dry clean. Do not iron.

Special maintenance

Special maintenance is required in case of failure or breakage following unforeseeable accidents or inappropriate use of the device.

The situations that might occur are entirely unpredictable; therefore, it is not possible to describe appropriate countermeasures.

If necessary, consult with the Moveo S.r.l. technical service for appropriate instructions for the specific situation.

All special maintenance activities must in any case be carried out by specialized and authorized personnel, under penalty of voiding the warranty.

Disposal

Any reuse of parts of the medical device is the responsibility of the user. The materials with which the medical device is made do not require special disposal procedures. Refer to local regulations for the disposal of unsorted residual waste. Remove the plastic and metal parts from the medical device and dispose of them as prescribed by law.

Do not dispose of or abandon the product in the environment for any reason whatsoever.

NOTE

THE MANUFACTURER IS NOT LIABLE FOR DAMAGE CAUSED BY THE DEVICE IF THIS IS NOT USED IN ITS COMPLETE CONFIGURATION AND FOR ITS INTENDED PURPOSE AS SPECIFIED IN THIS MANUAL.

THE MANUFACTURER IS NOT LIABLE FOR ANY PERSONAL INJURY OR DAMAGE TO PROPERTY RESULTING FROM THE RECOVERY AND REUSE OF PARTS ONCE THE DEVICE HAS BEEN DISMANTLED.



Declaration of conformity

The undersigned Mr. Fausto Antonio Panizzolo

in his capacity as the legal representative of Moveo S.r.l.

with registered office in: Via Monsignor Fortin 37/38, 35128 Padua, Italy

VAT number: 05236760285

declares

that the product: hip orthosis

Model and code: ExoBand

Medical device class: I

Date of manufacture: January 2022

Lot No: _____

was manufactured in compliance with the following directives and standards:

- Regulation (EU) no. 745/2017 on medical devices, amending Directive 2001/83/EC, Regulation (EC) no. 178/2002 and Regulation (EC) no. 1223/2009 and repealing Council Directives 90/385/EEC and 93/42/EEC
- IEC 61882:2016 Standard: Hazard and operability studies (HAZOP studies) - Application guide
- EN 980:2009 Symbols for use in the labelling of medical devices
- EN ISO 15223-1: 2012 Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements
- EN 1041:2009- Information supplied by the manufacturer of medical devices
- EN 14971:2012 - Application of risk management to medical devices
- Directive 93/42/EEC, Leg. Decree 24/02/97 no. 46, Leg. Decree 25/02/98 no. 95 Council Directive on medical devices - Implementation of Directive 93/42/EEC on medical devices. Changes to Leg. Decree 24 February 1997 no. 46 (which will coexist with the new Regulation until 2020)
- Directive 2007/47/EEC Leg. Decree 25/01/10 no. 37 Council Directive on medical devices - Implementation of Directive 93/42/EEC on medical devices. (which will coexist with the new Regulation until 2020)
- Directive 2001/95/EC also known as the general product safety directive
- Regulation (EU) No 1007/2011 of the European Parliament and of the Council of 27 September 2011 on textile fibre names and related labelling and marking of the fiber composition of textile products.

It therefore complies with current directives and standards.

This declaration of conformity is issued under the sole responsibility of the manufacturer.

Place: Padova

Date: 23 December 2019

Rev. 1

Signature:



Via Lauro, 95 Cadoneghe Padova Italy

Consulenti e Periti per: Guardia di Finanza, Tributaria, Autorità Doganali, Carabinieri, Polizia di Stato, Unioncamere, Courts of law.

The proper composition of the technical file as well as the documents prepared by C&C s.a.s. have been verified by Renato Carraro, engineer.

Basic UDI-DI: 805934045EXOBAND0003N